

Melamine Detection Test Service from Vet Food Agro Diagnostics (M) Sdn Bhd

Our laboratory is offering Melamine Detection Test and the details of the test are as follows:

Type of test: ELISA for Melamine Detection

Matrices: gluten, pet food (wet and dry)

Sample size: Send us 100g of each sample for testing (minimum of at least 50g)

Range of quantitation: Gluten: 10ppm to 250 ppm
Pet food (wet): 2ppm to 50 ppm
Pet food (dry): 4ppm to 100 ppm

Do call us for further enquiries at 03 – 7873 6405