



Highly Sensitive Melamine Test with Liquid Chromatography Tandem Mass Spectrometry (LC/MSMS)

Vet Food Agro Diagnostics (M) Sdn Bhd is proud to offer our highly sensitive and specific melamine test for animal feed, yeast, fish meal, soy bean meal, and milk powder.

Our melamine test is based on LC/MSMS method that produces reliable quantitative result with the Limit of Detection (LOD) down to 0.025 ppm which is 10 times lower than the current Maximum Residual Limit (MRL) of the Ministry of Health (MOH) at the level of 0.25 ppm. With such high sensitivity you can be rest assured that your raw materials/ products are free from melamine contamination.

For all your melamine testing requirements, please do not hesitate to contact us at:

Vet Food Agro Diagnostics (M) Sdn Bhd (560721-X)
(SAMM TESTING Laboratory No. 374)
Lot 18B, Jalan 241, Sekysen 51A,
46100 Petaling Jaya, Selangor, Malaysia.
P.O Box 150, 46710 Petaling Jaya, Selangor, Malaysia.
Tel: 603-7873 6405/ 7873 7355 (Ext 358)
Fax: 603-7873 9209